



***Ensure Patient Safety with
High Level Instrument
Disinfection***

Glutramil[®] OPA is a fast & effective high-level disinfectant for reprocessing reusable heat-sensitive devices such as endoscopes, respiratory therapy and anesthesia equipments.

- Excellent broad spectrum microbicidal activity: bactericidal, fungicidal, virucidal, sporicidal, tuberculocidal.
- Fast turn around time; provides high-level disinfection in 12 min.
- Easy to Use - No activation or mixing needed.
- Can be used/reused up to 14 days.
- Non corrosive solution.
- Easily rinses from equipment, especially lumened devices.
- Excellent material compatibility and can therefore be used to disinfect a wide range of medical instruments made of aluminum, brass, copper, stainless steel, plastics, elastomers, and dental materials.

COMPOSITION:

Each 100gm contains	
Ortho-phthalaldehyde	0.55 % w/w
Purified Water IP	q.s.

HOW TO USE :

Immerse clean, dry instruments completely in **Glutramil[®] OPA**. Soak the instruments for 12 min at 20°C or when using automatic endoscope reprocessor, for 5 minutes at 25°C to destroy pathogenic micro-organisms to achieve high level disinfection. Following disinfection, rinse instruments thoroughly, flushing the channels with portable/sterile water. Ensure thorough rinsing to remove any residual disinfectant.

PACK SIZE:

1 Litre
5 Litre

